



STERILE ZYCEINE 4YPINS (3-PART POUCH)

HYPOCHLOROUS ACID

ITEM NO. 0139 (5 X 25 WIPES PER CASE)

EFFECTIVE CLEANROOM SPORICIDES ARE USUALLY AGGRESSIVE CHEMICALS AND WHILST THEY CAN GIVE GOOD MICROBIAL EFFICACY THIS CAN OFTEN LEAD TO CORROSION PROBLEMS AND DAMAGE TO EXPENSIVE EQUIPMENT.

ZYCEINE 4YPINS IS A UNIQUE SPORICIDE IN A 3-PART POUCH THAT WHEN ACTIVATED GIVES GOOD MICROBIAL EFFICACY. ZYCEINE IS A TOTALLY STABLE PRODUCT UNTIL ACTIVATED.

ZYCEINE HAS A PROVEN FUNGICIDAL EFFICACY OF GREATER THAN A 3-LOG REDUCTION IN 2 MINUTES, AND A YEASTICIDAL EFFICACY GREATER THAN A 4-LOG REDUCTION IN 2 MINUTES.

- FAST ACTING SPORICIDE.
- EFFICACY INDEPENDENTLY VALIDATED.
- MADE WITH WFI-QUALITY WATER.
- 0.2-MICRON FILTERED.
- MANUFACTURED IN GRADE C/D CLEANROOMS.
- RE-SEALABLE TO PREVENT DRYING OUT.
- 30 GSM POLYPROPYLENE 25CM X 25CM WIPES.
- BROAD SPECTRUM ACTIVITY.
- FULLY TRACEABILITY WITH SUPPORTING ANALYSIS.
- GAMMA IRRADIATED FOR STERILITY.
- 2-YEAR UNOPENED SHELF LIFE.
- ALSO AVAILABLE IN 200ML/1L FLUID & 5 WIPES IN 2-PART POUCH FORMATS.



DIRECTION FOR USE:

TO ACTIVATE: 1. HOLD AND ROLL THE FLUID SECTION (1ST POUCH) INWARDS TOWARDS THE SEAL, PRESSING UNTIL THE VULNERABLE SEAL BREAKS, ALLOWING THE FLUID TO ENTER THE MIDDLE POUCH. 2. MANIPULATE CONTENTS VIGOROUSLY TO HOMOGENISE FOR 30 SECONDS. 3. LEAVE FOR 5 MINUTES AND AGITATE BRIEFLY BEFORE ROLLING THE MIXED SOLUTION TOWARDS THE LARGER SECTION AND BREAKING THE SECOND SEAL. 4. ONCE THE FLUID IS TRANSFERRED TO THE WIPES, AGITATE VIGOROUSLY BY PRESSING THE WIPES IN ALL DIRECTIONS FOR 30 SECONDS AND THEN LEAVE TO REST FOR 5 MIN. 5. OPEN AND RECLOSE THE POUCH AT THE ZIP SEAL. 6. USE WITHIN 8 HOURS OF ACTIVATION.

Use Biocides Safely. Always read the label and product information.

THE INFORMATION CONTAINED IN THIS DATA SHEET IS ACCURATE TO THE BEST OF OUR KNOWLEDGE AND BELIEF. NO WARRANTY IS, HOWEVER, GIVEN OR IMPLIED. CONSULT AGMA LTD CONCERNING ALL NON-STANDARD APPLICATIONS.

ISSUE NO. 1-02/26