



STERILE ALCOHOL 900ML TRIGGER SPRAY DE 70% IN WFI, BAG IN BOTTLE TRIGGER SPRAY. ITEM NO. 105665

STERILE DE 70% IN WFI IS A PROVEN MEANS FOR CONTINUAL DISINFECTION WITHIN PHARMACEUTICAL & HOSPITAL CLEANROOMS.

BAG IN BOTTLE TRIGGER SPRAYS ARE COMMONLY USED IN CLEANROOMS BOTH IN HOSPITALS AND PHARMACEUTICAL MANUFACTURING FACILITIES. THEY ARE A HIGH COST-EFFECTIVE MEANS OF ROUTINE DISINFECTION WHEN USED WITH STERILE DRY WIPES.

AS WELL AS A PROVEN DISINFECTANT, DE 70% IN WFI EVAPORATES FULLY LEAVING NO RESIDUES.

- EFFICACY INDEPENDENTLY VALIDATED.
- WFI QUALITY WATER – HIGH PURITY WATER.
- ERGONOMIC TRIGGER WITH EFFECTIVE RETURN MECHANISM.
- BAG IN BOTTLE TO MAINTAIN STERILITY IN USE.
- 13 WEEK IN-USE VALIDATION.
- MANUFACTURED IN ISO CLASS 5 CLEANROOM.
- FULLY TRACEABILITY WITH SUPPORTING ANALYSIS.
- GAMMA IRRADIATED FOR STERILITY.
- 2-YEAR SHELF LIFE.
- DOUBLE BAGGED.
- AVAILABLE IN STERILE WET WIPES, AEROSOLS AND 5L FORMATS.



DIRECTION FOR USE:

AGMA DE 70% IN WFI IS MANUFACTURED FOR USE ON CLEANROOM CEILINGS, WALLS, FLOORS AND EQUIPMENT HARD SURFACES.

ENSURE COMPLETE WETTING OF THE SURFACE. DO NOT USE TO TREAT SURFACES THAT ARE PARTICULARLY SENSITIVE TO ALCOHOL SUCH AS ACRYLIC GLASS. ALWAYS CLOSE CONTAINER AFTER USE.

Use Biocides Safely. Always read the label and product information.

THE INFORMATION CONTAINED IN THIS DATA SHEET IS ACCURATE TO THE BEST OF OUR KNOWLEDGE AND BELIEF. NO WARRANTY IS, HOWEVER, GIVEN OR IMPLIED. CONSULT AGMA LTD CONCERNING ALL NON-STANDARD APPLICATIONS.

ISSUE NO. 2-10/22