



STERILE ALCOHOL 900ML TRIGGER SPRAY

IPA 70% WFI BAG IN BOTTLE TRIGGER SPRAY

ITEM NO. 105666 (6 X 900ML)

STERILE IPA 70% WFI ARE A PROVEN MEANS FOR CONTINUAL DISINFECTION WITHIN PHARMACEUTICAL & HOSPITAL CLEANROOMS.

BAG IN BOTTLE TRIGGER SPRAYS ARE COMMONLY USED IN CLEANROOMS BOTH IN HOSPITALS AND PHARMACEUTICAL MANUFACTURING FACILITIES. THEY ARE A HIGH COST-EFFECTIVE MEANS OF ROUTINE DISINFECTION WHEN USED WITH STERILE DRY WIPES.

AS WELL AS A PROVEN DISINFECTANT IPA 70% EVAPORATES FULLY LEAVING NO RESIDUES.

- EFFICACY INDEPENDENTLY VALIDATED.
- WFI QUALITY WATER – HIGH PURITY WATER.
- 0.2 MICRON FILTERED.
- ERGONOMIC TRIGGER WITH EFFECTIVE RETURN MECHANISM.
- BAG IN BOTTLE TO MAINTAIN STERILITY IN USE.
- 13 WEEK IN-USE VALIDATION.
- MANUFACTURED IN GRADE C/D CLEANROOMS AND FILLED UNDER GRADE A CONDITIONS.
- FULLY TRACEABILITY WITH SUPPORTING ANALYSIS.
- GAMMA IRRADIATED FOR STERILITY.
- 2 YEAR SHELF LIFE.
- DOUBLE BAGGED.
- AVAILABLE WITH ALTERNATE ADJUSTABLE TRIGGER HEAD, 5L CAPPED BOTTLES AND WIPE FORMATS.



DIRECTION FOR USE:

AGMA IPA 70% IS MANUFACTURED TO HELP CLEANROOM OPERATORS MEET THE STANDARDS OF ANNEX 1, FOR USE ON CLEANROOM CEILINGS, WALLS, FLOORS, AND EQUIPMENT HARD SURFACES. ENSURE COMPLETE WETTING OF THE SURFACE. DO NOT USE TO TREAT SURFACES THAT ARE PARTICULARLY SENSITIVE TO ALCOHOL SUCH AS ACRYLIC GLASS. ALWAYS CLOSE THE CONTAINER AFTER USE.

Use Biocides Safely. Always read the label and product information.

THE INFORMATION CONTAINED IN THIS DATA SHEET IS ACCURATE TO THE BEST OF OUR KNOWLEDGE AND BELIEF. NO WARRANTY IS, HOWEVER, GIVEN OR IMPLIED. CONSULT AGMA LTD. CONCERNING ALL NON-STANDARD APPLICATIONS.

ISSUE NO.1-02/25