



# DE 70% 5L CAPPED BOTTLE

DE 70% IN WFI

ITEM NO. 1116 (4 X 5L PER CASE)

0.2 MICRON FILTERED

AGMA DE 70% IN WFI IS A PROVEN MEANS FOR CONTINUAL DISINFECTION WITHIN CLEANROOMS.

AGMA DE 70% IN WFI IS COMMONLY USED IN CLEANROOM FACILITIES INCLUDING IN HOSPITAL OR PHARMACEUTICAL ENVIRONMENTS. AGMA DE 70% IN WFI IS A RELIABLE, COST-EFFECTIVE MEANS OF ROUTINE DISINFECTION WHEN USED WITH DRY MOPS OR WIPES.

AS WELL AS A PROVEN DISINFECTANT, DE 70% EVAPORATES FULLY LEAVING NO RESIDUES.

- EFFICACY INDEPENDENTLY VALIDATED.
- 0.2 MICRON FILTERED.
- MADE WITH WFI QUALITY WATER.
- MANUFACTURED IN GRADE C/D CLEANROOMS.
- IDEAL FOR GRADE C/D CLEANROOMS WHERE A CONTAMINATION CONTROL STRATEGY HAS BEEN UNDERTAKEN.
- FULLY TRACEABILITY WITH SUPPORTING ANALYSIS.
- 2 YEAR SHELF LIFE.
- AVAILABLE IN TRIGGER SPRAY & WET WIPE FORMATS.



## DIRECTION FOR USE:

AGMA DE 70% IS MANUFACTURED TO HELP CLEANROOM OPERATORS MEET THE STANDARDS OF ANNEX 1, FOR USE ON CLEANROOM CEILINGS, WALLS, FLOORS AND EQUIPMENT HARD SURFACES. ENSURE COMPLETE WETTING OF THE SURFACE. DO NOT USE TO TREAT SURFACES THAT ARE PARTICULARLY SENSITIVE TO ALCOHOL SUCH AS ACRYLIC GLASS. ALWAYS CLOSE CONTAINER AFTER USE.

Use Biocides Safely. Always read the label and product information.

THE INFORMATION CONTAINED IN THIS DATA SHEET IS ACCURATE TO THE BEST OF OUR KNOWLEDGE AND BELIEF. NO WARRANTY IS, HOWEVER, GIVEN OR IMPLIED. CONSULT AGMA LTD CONCERNING ALL NON-STANDARD APPLICATIONS.

ISSUE NO. 1-03/23