



IPA 70% IN WFI 5L BOTTLE

ITEM NO. 1214 (4 X 5L CAP BOTTLES PER CASE)

0.2 MICRON FILTERED

AGMA IPA 70% IN WFI IS A PROVEN MEANS FOR CONTINUAL DISINFECTION WITHIN CLEANROOMS.

AGMA IPA 70% IN WFI IS COMMONLY USED IN CLEANROOM FACILITIES INCLUDING IN HOSPITAL OR PHARMACEUTICAL ENVIRONMENTS. AGMA IPA 70% IN WFI IS A RELIABLE, COST-EFFECTIVE MEANS OF ROUTINE DISINFECTION WHEN USED WITH DRY MOPS OR WIPES.

AS WELL AS A PROVEN DISINFECTANT, IPA 70% EVAPORATES FULLY LEAVING NO RESIDUES.

- MANUFACTURED IN GRADE C/D CLEANROOM
- IDEAL FOR GRADE C/D CLEANROOMS WHERE A CONTAMINATION CONTROL STRATEGY HAS BEEN UNDERTAKEN.
- 0.2 MICRON FILTERED
- HIGH PURITY WFI QUALITY WATER
- FULLY TRACEABILITY WITH SUPPORTING ANALYSIS.
- 2 YEAR SHELF LIFE.
- EFFICACY VALIDATED EXTERNALLY
- AVAILABLE IN TRIGGER SPRAY & WET WIPE FORMATS



DIRECTION FOR USE:

AGMA IPA 70% IS MANUFACTURED TO HELP CLEANROOM OPERATORS MEET THE STANDARDS OF ANNEX 1, FOR USE ON CLEANROOM CEILINGS, WALLS, FLOORS, AND EQUIPMENT HARD SURFACES. ENSURE COMPLETE WETTING OF THE SURFACE. DO NOT USE TO TREAT SURFACES THAT ARE PARTICULARLY SENSITIVE TO ALCOHOL SUCH AS ACRYLIC GLASS. ALWAYS CLOSE THE CONTAINER AFTER USE.

Use Biocides Safely. Always read the label and product information.

THE INFORMATION CONTAINED IN THIS DATA SHEET IS ACCURATE TO THE BEST OF OUR KNOWLEDGE AND BELIEF. NO WARRANTY IS, HOWEVER, GIVEN OR IMPLIED. CONSULT AGMA LTD CONCERNING ALL NON-STANDARD APPLICATIONS.

ISSUE NO. 1-03/25