



# STERILE IPA 70% IN WFI WIPES (DB)

68GSM PRE-IMPREGNATED – READY TO USE / EASY TO USE.  
ITEM NO. 0122 (10 X 10 WIPES)

STERILE IPA 70% WFI ARE A PROVEN MEANS FOR CONTINUAL DISINFECTION WITHIN PHARMACEUTICAL & HOSPITAL CLEANROOMS.

READY TO USE WIPES ARE A SIMPLE MEANS OF DISINFECTING A VARIETY OF SURFACES AND CLEANROOM ITEMS. THEY HAVE THE ADDED BENEFIT OF BEING ABLE TO PHYSICALLY REMOVE PARTICLES AND SPORES IN ADDITION TO PROVIDING GOOD DISINFECTION AGAINST VEGETATIVE MICROBES.

COST IN USE, TAKING INTO ACCOUNT OPERATOR TIME, GIVES FINANCIAL BENEFITS WHEN COMPARED TO DRY WIPES AND SPRAYS.

- READY TO USE-SAVES TIME.
- RE-SEALABLE TO PREVENT DRYING OUT.
- 0.2 MICRON FILTERED.
- 28CM X 38CM WIPES.
- 68GSM PER M<sup>2</sup> LOW PARTICULATE POLYCELLULOSE WIPE.
- BROAD SPECTRUM ACTIVITY.
- WFI QUALITY WATER – HIGH PURITY WATER.
- MANUFACTURED IN GRADE C/D CLEANROOMS.
- FULLY TRACEABILITY WITH SUPPORTING ANALYSIS.
- GAMMA IRRADIATED FOR STERILITY.
- DOUBLE BAGGED.
- 2 YEAR SHELF LIFE.



## DIRECTION FOR USE:

AGMA IPA 70% IS MANUFACTURED TO HELP CLEANROOM OPERATORS MEET THE STANDARDS OF ANNEX 1, FOR USE ON CLEANROOM CEILINGS, WALLS, FLOORS, AND EQUIPMENT HARD SURFACES. ENSURE COMPLETE WETTING OF THE SURFACE. DO NOT USE TO TREAT SURFACES THAT ARE PARTICULARLY SENSITIVE TO ALCOHOL SUCH AS ACRYLIC GLASS. ALWAYS CLOSE THE CONTAINER AFTER USE.

Use Biocides Safely. Always read the label and product information.

THE INFORMATION CONTAINED IN THIS DATA SHEET IS ACCURATE TO THE BEST OF OUR KNOWLEDGE AND BELIEF. NO WARRANTY IS, HOWEVER, GIVEN OR IMPLIED. CONSULT AGMA LTD. CONCERNING ALL NON-STANDARD APPLICATIONS.

ISSUE NO. 1-03/25