



## STERILE NEUTRAL DETERGENT

ITEM NO 192900: 10 X POUCHES OF 10 WIPES PER CASE

FOR EFFECTIVE DISINFECTION TO BE ACHIEVED IT IS IMPERATIVE THAT SURFACES ARE FREE FROM ANY VISIBLE SOILING OR BUILD-UP OF RESIDUE FROM PREVIOUS DISINFECTANT APPLICATIONS.

AGMA HAVE CONSIDERABLE EXPERIENCE IN SURFACTANT TECHNOLOGY AND STERILE CLEANING / DISINFECTION. WE HAVE COMBINED THIS EXPERTISE TO DEVELOP A STERILE NEUTRAL DETERGENT SPECIFICALLY FOR USE IN CLEANROOMS AND PHARMACIES.

WITH HIGHLY EFFECTIVE SURFACE CLEANING PROPERTIES, AGMA STERILE NEUTRAL DETERGENT IS IDEAL WHEN A DEEP CLEAN IS REQUIRED. FORMULATED TO BE LOW FOAMING IT DOES NOT LEAVE SIGNIFICANT RESIDUES OR SMEARS. IN ADDITION, ITS PH HELPS PROTECT MORE DELICATE WORK SURFACES AND THOSE PRONE TO CORROSION.

- MANUFACTURED IN CLEANROOMS.
- MADE WITH WFI QUALITY WATER.
- 0.2 MICRON FILTERED.
- GAMMA IRRADIATED FOR STERILITY
- LOW FOAMING.
- HIGHLY EFFECTIVE SURFACE CLEANING PROPERTIES.
- REMOVES GREASES, DISINFECTANT RESIDUE ETC.
- DOES NOT LEAVE SIGNIFICANT RESIDUES.
- 69GSM, LOW LINTING POLYCELLULOSE WIPE
- 28CM X 38CM FOR LARGE SURFACE AREA COVERAGE
- FULLY TRACEABLE WITH SUPPORTING ANALYSIS.
- 2 YEAR SHELF LIFE.



### DIRECTION FOR USE:

AGMA NEUTRAL DETERGENT IS MANUFACTURED FOR USE ON CLEANROOM HARD SURFACES TO REMOVE PARTICULATE CONTAMINATION SUCH AS DISINFECTANT RESIDUE. USE THESE WIPES TO PERFORM AN INITIAL DEEP CLEAN PRIOR TO THE USE OF ANY DISINFECTANT.

THE INFORMATION CONTAINED IN THIS DATA SHEET IS ACCURATE TO THE BEST OF OUR KNOWLEDGE AND BELIEF. NO WARRANTY IS, HOWEVER, GIVEN OR IMPLIED. CONSULT AGMA LTD CONCERNING ALL NON-STANDARD APPLICATIONS.

ISSUE NO. 4-09/23