



## STERILE NEUTRAL DETERGENT RFU 5 L

2 X 5L CAPPED BOTTLES PER CASE

ITEM NO 0008

FOR EFFECTIVE DISINFECTION TO BE ACHIEVED IT IS IMPERATIVE THAT SURFACES ARE FREE FROM ANY VISIBLE SOILING OR BUILD-UP OF RESIDUE FROM PREVIOUS DISINFECTANT APPLICATIONS.

AGMA HAVE CONSIDERABLE EXPERIENCE IN SURFACTANT TECHNOLOGY AND CLEANING / DISINFECTION. WE HAVE COMBINED THIS EXPERTISE TO DEVELOP A NEUTRAL DETERGENT SPECIFICALLY FOR USE IN CLEANROOMS AND PHARMACIES.

WITH HIGHLY EFFECTIVE SURFACE CLEANING PROPERTIES IT IS IDEAL WHEN A DEEP CLEAN IS REQUIRED. FORMULATED TO BE LOW FOAMING IT DOES NOT LEAVE SIGNIFICANT RESIDUES OR SMEARS. IN ADDITION, IT IS PH NEUTRAL BEFORE IRRADIATION, SO HELPS TO PROTECT MORE DELICATE WORK SURFACES AND THOSE PRONE TO CORROSION.

- MADE WITH WFI QUALITY WATER.
- GAMMA IRRADIATED FOR STERILITY.
- 0.2 MICRON FILTERED.
- LOW FOAMING.
- HIGHLY EFFECTIVE SURFACE CLEANING PROPERTIES.
- REMOVES GREASES, DISINFECTANT RESIDUE.
- MANUFACTURED IN GRADE C/D CLEANROOMS.
- DOES NOT LEAVE SIGNIFICANT RESIDUES.
- FULLY TRACEABLE WITH SUPPORTING ANALYSIS.
- 2 YEAR SHELF LIFE.
- AVAILABLE IN VARIOUS OTHER FORMATS.



### DIRECTION FOR USE

AGMA NEUTRAL DETERGENT IS MANUFACTURED TO HELP CLEANROOM OPERATORS MEET THE STANDARDS OF ANNEX 1. IDEAL FOR USE ON CLEANROOM HARD SURFACES TO REMOVE RESIDUES AND FOR USE BETWEEN ROTATING OF DISINFECTANTS. AN INITIAL DEEP CLEAN SHOULD BE PERFORMED PRIOR TO THE USE OF ANY DISINFECTANT TO HELP THE REMOVAL OF ANY SOILS OR HEAVY RESIDUE BUILD UP.

THE INFORMATION CONTAINED IN THIS DATA SHEET IS ACCURATE TO THE BEST OF OUR KNOWLEDGE AND BELIEF. NO WARRANTY IS, HOWEVER, GIVEN OR IMPLIED. CONSULT AGMA LTD CONCERNING ALL NON-STANDARD APPLICATIONS.

ISSUE NO. 1-02/25