

PRODUCTS

Sterile Alcohol Products

Alcohol-based products are typically available as both 70% Isopropanol (IPA) & 70% Denatured Ethanol (DE). All AGMA alcohol products use European Pharmacopeia-grade alcohol before being blended with injection-quality water in order to guarantee endotoxin levels below 0.25EU/ml, prior to being sterilised by a validated gamma irradiation process at no less than 25kGy, (5L product 20kGy - 45kGy). All products are available with the necessary supporting documentation required to demonstrate a premium quality product that is suitable for use in all cleanroom areas as required. All CEN test data is available on request.


IPA 70%: EN14476 | EN1276 | EN1650 | EN13697 Validated

DE 70%: EN14476 | EN1276 | EN1650 | EN13697 | EN1275 | EN1040 Validated



IPA & DE 70% 5L

- Ø IPA & DE 70% in WFI, guaranteeing a low endotoxin level and suitability for product contact areas.
- Ø 0.2 micron filtered and 20kGy - 45kGy gamma irradiated for sterility.
- Ø Batch coded with Detex dot to clearly indicate successful irradiation.
- Ø 2 year unopened shelf life.



PRODUCT INFO

Item No.	Description	Content	Packaging Format
0036	Sterile DE 70% in WFI	5L	2 x 5L individually double bagged
0038	Sterile IPA 70% in WFI	5L	2 x 5L individually double bagged