



## STERILE NEUTRAL DETERGENT

## ITEM NO 0055: 20 X 500ML CAPPED BOTTLES PER CASE

FOR EFFECTIVE DISINFECTION TO BE ACHIEVED IT IS IMPERATIVE THAT SURFACES ARE FREE FROM ANY VISIBLE SOILING OR BUILD-UP OF RESIDUE FROM PREVIOUS DISINFECTANT APPLICATIONS.

AGMA HAVE CONSIDERABLE EXPERIENCE IN SURFACTANT TECHNOLOGY AND STERILE CLEANING / DISINFECTION. WE HAVE COMBINED THIS EXPERTISE TO DEVELOP A STERILE NEUTRAL DETERGENT SPECIFICALLY FOR USE IN CLEANROOMS AND PHARMACIES.

WITH HIGHLY EFFECTIVE SURFACE CLEANING PROPERTIES, AGMA STERILE NEUTRAL DETERGENT IS IDEAL WHEN A DEEP CLEAN IS REQUIRED. FORMULATED TO BE LOW FOAMING IT DOES NOT LEAVE SIGNIFICANT RESIDUES OR SMEARS. IN ADDITION, ITS PH HELPS PROTECT MORE DELICATE WORK SURFACES AND THOSE PRONE TO CORROSION.

- MANUFACTURED IN CLEANROOMS.
- MADE WITH WFI OUALITY WATER.
- 0.2 MICRON FILTERED.
- INDIVIDUALLY DOUBLE BAGGED
- GAMMA IRRADIATED FOR STERILITY
- LOW FOAMING.
- HIGHLY EFFECTIVE SURFACE CLEANING PROPERTIES.
- REMOVES GREASES, DISINFECTANT RESIDUE.
- DOES NOT LEAVE SIGNIFICANT RESIDUES.
- FULLY TRACEABLE WITH SUPPORTING ANALYSIS.
- 2 YEAR SHELF LIFE.
- AVAILABLE IN VARIOUS OTHER FORMATS



## **DIRECTION FOR USE:**

AGMA NEUTRAL DETERGENT IS MANUFACTURED FOR USE ON CLEANROOM HARD SURFACES TO REMOVE PARTICULATE CONTAMINATION SUCH AS DISINFECTANT RESIDUE. USE WITH WIPES TO PERFORM AN INITIAL DEEP CLEAN PRIOR TO THE USE OF ANY DISINFECTANT. BEING MANUFACTURED TO ENABLE COMPLIANCE WITH ANNEX 1, ITEM 0055 IS PERFECT FOR USE IN GRADE A CLEANROOMS OR LOWER. USE WITH WIPES TO PERFORM AN INITIAL DEEP CLEAN PRIOR TO THE USE OF ANY DISINFECTANT

THE INFORMATION CONTAINED IN THIS DATA SHEET IS ACCURATE TO THE BEST OF OUR KNOWLEDGE AND BELIEF. NO WARRANTY IS, HOWEVER, GIVEN OR IMPLIED.

CONSULT AGMA LTD CONCERNING ALL NON-STANDARD APPLICATIONS.

ISSUE NO. 4-09/23

