



DE 70% 500 ML TRIGGER SPRAY

DE 70% IN WFI

ITEM NO. 1139

0.2 MICRON FILTERED

DE 70% WFI IS A PROVEN MEANS FOR CONTINUAL DISINFECTION WITHIN PHARMACEUTICAL AND HOSPITAL CLEANROOMS.

DE 70% ALCOHOL TRIGGER SPRAYS ARE A HIGHLY COST-EFFECTIVE MEANS OF ROUTINE DISINFECTION WHEN USED WITH STERILE DRY WIPES.

AS WELL AS A PROVEN DISINFECTANT DE 70% EVAPORATES FULLY SO THERE IS NO RESIDUE.

- EFFICACY INDEPENDENTLY VALIDATED.
- MADE WITH WFI QUALITY WATER.
- ERGONOMIC TRIGGER WITH EFFECTIVE RETURN MECHANISM.
- 0.2 MICRON FILTERED.
- RELIABLE BOTTLE WITH HEADROOM FOR EXPANSION.
- MANUFACTURED IN GRADE C/D CLEANROOMS.
- IDEAL FOR GRADE C/D CLEANROOMS WHERE A CONTAMINATION CONTROL STRATEGY HAS BEEN UNDERTAKEN.
- FULL TRACEABILITY WITH SUPPORTING DOCUMENTATION.
- CONSISTENT SPRAY PATTERN.
- ALSO AVAILABLE IN 1L TRIGGER SPRAY.



DIRECTION FOR USE:

AGMA DE 70% IS MANUFACTURED TO HELP CLEANROOM OPERATORS MEET THE STANDARDS OF ANNEX 1, FOR USE ON CLEANROOM CEILINGS, WALLS, FLOORS AND EQUIPMENT HARD SURFACES. ENSURE COMPLETE WETTING OF THE SURFACE. DO NOT USE TO TREAT SURFACES THAT ARE PARTICULARLY SENSITIVE TO ALCOHOL SUCH AS ACRYLIC GLASS. ALWAYS CLOSE CONTAINER AFTER USE.

Use Biocides Safely. Always read the label and product information.

THE INFORMATION CONTAINED IN THIS DATA SHEET IS ACCURATE TO THE BEST OF OUR KNOWLEDGE AND BELIEF. NO WARRANTY IS, HOWEVER, GIVEN OR IMPLIED. CONSULT AGMA LTD CONCERNING ALL NON-STANDARD APPLICATIONS. ISSUE NO. 1-03/25