



IPA 70% 1L TRIGGER SPRAY

IPA 70% IN WFI

ITEM NO. 1215 (15 X 1L)

0.2 MICRON FILTERED

IPA 70% IN WFI IS A PROVEN MEANS FOR CONTINUAL DISINFECTION WITHIN PHARMACEUTICAL & HOSPITAL CLEANROOMS.

AGMA IPA 70% TRIGGER SPRAYS ARE A HIGH QUALITY, COST- EFFECTIVE MEANS OF ROUTINE DISINFECTION, ESPECIALLY WHEN USED WITH DRY WIPES.

AS WELL AS A PROVEN DISINFECTANT IPA 70% EVAPORATES FULLY SO THERE IS NO RESIDUE.

- MANUFACTURED IN GRADE C/D CLEANROOMS.
- IDEAL FOR GRADE C/D CLEANROOMS WHERE A CONTAMINATION CONTROL STRATEGY HAS BEEN UNDERTAKEN.
- EFFICACY INDEPENDENTLY VALIDATED.
- 0.2 MICRON FILTERED
- HIGH PURITY, WFI QUALITY WATER
- ERGONOMIC TRIGGER WITH EFFECTIVE RETURN MECHANISM.
- RELIABLE BOTTLE.
- FULL TRACEABILITY WITH SUPPORTING ANALYSIS.
- CONSISTENT SPRAY PATTERN
- ALSO AVAILABLE IN 500ML TRIGGERS, 5L AND WET WIPE FORMATS



DIRECTION FOR USE:

AGMA IPA 70% IS MANUFACTURED TO HELP CLEANROOM OPERATORS MEET THE STANDARDS OF ANNEX 1, FOR USE ON CLEANROOM CEILINGS, WALLS, FLOORS, AND EQUIPMENT HARD SURFACES. ENSURE COMPLETE WETTING OF THE SURFACE. DO NOT USE TO TREAT SURFACES THAT ARE PARTICULARLY SENSITIVE TO ALCOHOL SUCH AS ACRYLIC GLASS. ALWAYS CLOSE THE CONTAINER AFTER USE.

Use Biocides Safely. Always read the label and product information.

THE INFORMATION CONTAINED IN THIS DATA SHEET IS ACCURATE TO THE BEST OF OUR KNOWLEDGE AND BELIEF. NO WARRANTY IS, HOWEVER, GIVEN OR IMPLIED. CONSULT AGMA LTD CONCERNING ALL NON-STANDARD APPLICATIONS.

ISSUE NO. 1-03/25