



# STERILE NEUTRAL DETERGENT MOP WIPES

## 20 X 5 PACKS PER CASE IN WFI PRE-IMPREGNATED – READY TO USE ITEM 0085

FOR EFFECTIVE DISINFECTION TO BE ACHIEVED IT IS IMPERATIVE THAT SURFACES ARE FREE FROM ANY VISIBLE SOILING OR BUILD-UP OF RESIDUE FROM PREVIOUS DISINFECTANTS.

READY TO USE STERILE NEUTRAL DETERGENT MOP WIPES ARE HIGHLY EFFECTIVE AT REMOVING GREASES, SPILLS AND DISINFECTANT RESIDUE FROM LARGER SURFACE AREAS.

OUR PRESATURATED SINGLE-USE MOP WIPES ARE READY TO USE AND PERFECT FOR A BUCKETLESS SYSTEM, SAVING OPERATOR TIME. THE MOP WIPES ARE EASY TO INSTALL AND DISPOSE OF AND ARE COMPATIBLE WITH WIDELY AVAILABLE 'PUSH CLIP' MOP HEADS.

- REMOVES GREASES, SPILLS AND DISINFECTANT RESIDUE.
- LOW FOAMING.
- NO SIGNIFICANT RESIDUES.
- FORMULATED TO BE PH NEUTRAL BEFORE IRRADIATION.
- MADE WITH WFI QUALITY WATER.
- MANUFACTURED IN GRADE C/D CLEANROOMS.
- POLYCELLULOSE 60GSM HIGH QUALITY, LOW LINTING MOP WIPE.
- 25CM X 50CM, ALLOWING FOR LARGE AREA COVERAGE.
- FULLY TRACEABLE WITH SUPPORTING ANALYSIS.
- GAMMA IRRADIATION FOR STERILITY.
- 20 MOP WIPES PER PACK, 5 PACKS PER CASE.
- 2 YEAR UNOPENED SHELF LIFE.



### DIRECTION FOR USE

AGMA NEUTRAL DETERGENT WIPES ARE MANUFACTURED TO HELP CLEANROOM OPERATORS MEET THE STANDARDS OF ANNEX 1. IDEAL FOR USE ON CLEANROOM HARD SURFACES TO REMOVE RESIDUES AND FOR USE BETWEEN ROTATING OF DISINFECTANTS. AN INITIAL DEEP CLEAN SHOULD BE PERFORMED PRIOR TO THE USE OF ANY DISINFECTANT TO HELP THE REMOVAL OF ANY SOILS OR HEAVY RESIDUE BUILD UP.

THE INFORMATION CONTAINED IN THIS DATA SHEET IS ACCURATE TO THE BEST OF OUR KNOWLEDGE AND BELIEF. NO WARRANTY IS, HOWEVER, GIVEN OR IMPLIED. CONSULT AGMA LTD CONCERNING ALL NON-STANDARD APPLICATIONS.

ISSUE NO. 1-02/25