



STERILE DE 70% 400ML AEROSOL

DE 70% IN WFI

ITEM NO. 105614 (16 X 400 ML)

STERILE DE 70% IN WFI 400ML AEROSOLS ARE AN EASY TO USE AND POPULAR MEANS OF APPLYING DISINFECTANT INSIDE CLEANROOMS AND FOR WIPING IN PROCEDURES.

THEIR LIGHT WEIGHT AND EASY TO OPERATE SPRAY ARE WELCOMED BY USERS WHO MAY FIND LARGER SPRAYS HEAVY AND REPETITIVE TRIGGER ACTION TIRING. ALSO GIVING A VERY REGULAR SPRAY PATTERN.

BEING PRESSURISED THERE IS NO ISSUE WITH SUCK - BACK AND CONTENTS WILL ALWAYS REMAIN STERILE IN USE.

- MADE WITH WFI QUALITY WATER.
- 0.2 MICRON FILTERED.
- MANUFACTURED IN GRADE C/D CLEANROOMS.
- LIGHTWEIGHT-VERY POPULAR WITH USERS.
- NO TIRING TRIGGER ACTION.
- PRESSURISED CONTAINER – ALWAYS MAINTAINS STERILITY.
- REGULAR SPRAY PATTERN.
- TRIGGER ATTACHMENT AVAILABLE.
- DOUBLE BAGGED.
- 3 YEAR SHELF LIFE.
- FULLY TRACEABLE WITH SUPPORTING ANALYSIS.



DIRECTION FOR USE:

AGMA DE 70% IS MANUFACTURED TO HELP CLEANROOM OPERATORS MEET THE STANDARDS OF ANNEX 1, FOR USE ON CLEANROOM CEILINGS, WALLS, FLOORS AND EQUIPMENT HARD SURFACES. ENSURE COMPLETE WETTING OF THE SURFACE. DO NOT USE TO TREAT SURFACES THAT ARE PARTICULARLY SENSITIVE TO ALCOHOL SUCH AS ACRYLIC GLASS. ALWAYS CLOSE CONTAINER AFTER USE.

Use Biocides Safely. Always read the label and product information.

THE INFORMATION CONTAINED IN THIS DATA SHEET IS ACCURATE TO THE BEST OF OUR KNOWLEDGE AND BELIEF. NO WARRANTY IS, HOWEVER, GIVEN OR IMPLIED. CONSULT AGMA LTD. CONCERNING ALL NON-STANDARD APPLICATIONS.

ISSUE NO. 1-02/25