



DE TUB WIPES

DE 70% V/V IN WFI
ITEM NO. 0140 (15 X 100 WIPES)

DENATURED ETHANOL (DE) 70% V/V IN WFI QUALITY WATER TUB WIPES ARE A PROVEN MEANS FOR CONTINUAL DISINFECTION WITHIN PHARMACEUTICAL & HOSPITAL CLEANROOMS.

READY TO USE FOR A VARIETY OF SURFACES, EQUIPMENT, AND DEVICES IN CLEANROOM AND CLINICAL ENVIRONMENTS. WIPES HAVE THE ADDED BENEFIT OF BEING ABLE TO PHYSICALLY REMOVE PARTICLES AND SPORES IN ADDITION TO PROVIDING GOOD DISINFECTION AGAINST VEGETATIVE MICROBES.

TAKING INTO ACCOUNT OPERATOR TIME, COST IN USE GIVES FINANCIAL BENEFITS WHEN COMPARED TO DRY WIPES AND SPRAYS.

- READY TO USE - SAVES TIME.
- RE-SEALABLE TO PREVENT DRYING OUT.
- 25GSM LOW PARTICULATE POLYPROPYLENE WIPE.
- 185 X 200 MM FOR LARGE SURFACE AREA COVERAGE
- BROAD SPECTRUM ACTIVITY.
- MADE WITH WFI QUALITY WATER.
- FULLY TRACEABILITY WITH SUPPORTING ANALYSIS.
- 0.2 MICRON FILTERED.
- 2 YEAR SHELF LIFE.
- DE IS AVAILALBE IN VARIOUS STERILE/ 0.2 FILTERED FORMATS.



DIRECTION FOR USE:

AGMA DE 70% V/V IS MANUFACTURED TO HELP CLEANROOM OPERATORS MEET THE STANDARDS OF ANNEX 1. USE FOR A VARIETY OF SURFACES, EQUIPMENT, AND DEVICES IN CLEANROOMS. DO NOT USE TO TREAT SURFACES THAT ARE PARTICULARLY SENSITIVE TO ALCOHOL SUCH AS ACRYLIC GLASS. ALWAYS CLOSE CONTAINER AFTER USE. ENSURE ADEQUATE VENTILATION.

Use Biocides Safely. Always read the label and product information

THE INFORMATION CONTAINED IN THIS DATA SHEET IS ACCURATE TO THE BEST OF OUR KNOWLEDGE AND BELIEF. NO WARRANTY IS, HOWEVER, GIVEN OR IMPLIED. CONSULT AGMA LTD. CONCERNING ALL NON-STANDARD APPLICATIONS.

ISSUE NO. 1-04/26