



STERILE DE 70% 400ML AEROSOL

DE 70% IN WFI

ITEM NO. 9711 (12 X 400 ML)

STERILE DE 70% IN WFI 400ML AEROSOLS ARE AN EASY TO USE AND POPULAR MEANS OF APPLYING DISINFECTANT INSIDE CLEANROOMS AND FOR WIPING IN PROCEDURES.

OFFERS AN ADVANCED BAG-ON-VALVE (BOV) AEROSOL DISPENSING SYSTEM WHERE A PRODUCT IS SEALED INSIDE A MULTI-LAYERED BAG, COMPLETELY SEPARATED FROM THE PROPELLANT. BOV TECHNOLOGY ALSO ALLOWS COMPLETE DECANTING OF FLUID.

360-DEGREE SPRAYER ALLOWS FOR A CONSISTENT SPRAY PATTERN AT ANY ANGLE. THIS REDUCES DOWNTIME, CHANCE OF WRIST FATIGUE AND ALLOWS AN EVEN APPLICATION OF THE PRODUCT.

- BAG-ON-VALVE TECHNOLOGY (BOV)
- COMPRESSED AIR USED AS PROPELLANT.
- MADE WITH WFI QUALITY WATER.
- 0.2 MICRON FILTERED.
- MANUFACTURED IN GRADE C/D CLEANROOMS.
- LIGHTWEIGHT – VERY POPULAR WITH USERS.
- PRESSURISED CONTAINER – ALWAYS MAINTAINS STERILITY.
- REGULAR SPRAY PATTERN.
- TRIGGER ATTACHMENT AVAILABLE.
- DOUBLE BAGGED.
- 3 YEAR SHELF LIFE.
- FULLY TRACEABLE WITH SUPPORTING ANALYSIS.



DIRECTION FOR USE:

AGMA DE 70% IN WFI QUALITY WATER IS MANUFACTURED TO HELP CLEANROOM OPERATORS MEET THE STANDARDS OF ANNEX 1, FOR USE ON CLEANROOM CEILINGS, WALLS, FLOORS AND EQUIPMENT HARD SURFACES. ENSURE COMPLETE WETTING OF THE SURFACE. DO NOT USE TO TREAT SURFACES THAT ARE PARTICULARLY SENSITIVE TO ALCOHOL SUCH AS ACRYLIC GLASS. ALWAYS CLOSE CONTAINER AFTER USE.

Use Biocides Safely. Always read the label and product information.

THE INFORMATION CONTAINED IN THIS DATA SHEET IS ACCURATE TO THE BEST OF OUR KNOWLEDGE AND BELIEF. NO WARRANTY IS, HOWEVER, GIVEN OR IMPLIED. CONSULT AGMA LTD. CONCERNING ALL NON-STANDARD APPLICATIONS.

ISSUE NO. 1-06/26