

PRODUCTS

Sterile Neutral Detergent

For effective disinfection to be achieved, it is imperative that surfaces are free from any visible soiling or build up of residue from previous disinfectant applications. AGMA have considerable expertise in surfactant technology and sterile cleaning/disinfection. We have combined this expertise to develop a Sterile Neutral Detergent specifically for use in cleanrooms. With highly effective surface cleaning properties, it is ideal when a deep clean is required. Formulated to be low foaming, it does not leave significant residues or smears. In addition, it is pH neutral before irradiation so helps protect more delicate work surfaces and those prone to corrosion.



STERILE NEUTRAL DETERGENT BOTTLES

- ∅ 500ml bottle for use over small areas, isolators etc.
- ∅ 1L & 5L, also available to suit requirements.
- ∅ Ideal for dosing sterile dry wipes and mops.
- ∅ 0.2 micron filtered and 25kGy gamma irradiated for sterility. 20-45kGy for 5L product.
- ∅ Manufactured using WFI quality water guaranteeing a low endotoxin level and suitability for product contact areas.
- ∅ Double bagged and batch coded with Detex dot to clearly indicate successful irradiation.
- ∅ 2 year unopened shelf life.



PRODUCT INFO

Item No.	Description	Content	Packaging Format
0055	Sterile Neutral Detergent in WFI Bottle	500ml	20 x 500ml capped bottle, individually double bagged
5572	Sterile Neutral Detergent in WFI Bottle	1L	15 x 1L, capped bottle, individually double bagged
0008	Sterile Neutral Detergent in WFI	5L	2 x 5L, capped bottle, individually double bagged