STERILE ALCOHOL 70% TRIGGER SPRAYS

- Trigger spray is optimised to help reduce excessive Ø strain on the operator.
- 'Bag in bottle' technology balances air pressure and Ø has validated in-use sterility for 13 weeks.
- 0.2 micron filtered and 25kGy gamma irradiated for Ø sterility.
- Available as IPA 70% in WFI or as DE 70% in WFI, Ø guaranteeing a low endotoxin level and suitability for product contact areas.
- Filled and double bagged under a Grade A Laminar Ø Air Flow in a Grade C (ISO Class 7) cleanroom.
- Batch coded with Detex dot to clearly indicate Ø successful irradiation.
- Ø 2 year unopened shelf life.

PRODUCT INFO

ltem No.	Description	Content	Packaging Format
105665	Sterile DE 70% in WFI Trigger Spray	900ml	6 x 900 ml Individually double bagged
105666	Sterile IPA 70% in WFI Trigger Spray	900ml	6 x 900 ml Individually double bagged



For validations see page 8



STERILE ALCOHOL 70% AEROSOL

Aerosols are an easy-to-use and popular means of applying sterile alcohol. Their light weight and easy-to-operate spray is welcomed by users who may find larger sprays heavy and the repetitive trigger action tiring. Being pressurised there is no issue with suck-back and contents will remain sterile in use. Available as IPA 70% in WFI or as DE 70% in WFI.

- Guaranteeing a low endotoxin level and suitability for product contact areas. Ø
- 0.2 micron filtered and 25kGy gamma irradiated for sterility. Ø
- Lightweight and easy to use. Ø
- No tiring trigger action. Ø
- Pressurised and airtight container maintains sterility. Ø
- Trigger attachments are available if required. Ø
- Double bagged for ease of transfer. Ø
- Batch coded with Detex dot to clearly indicate Ø successful irradiation.
- Incorporates a non-flammable, inert propellant. Ø
- Ø 3 year unopened shelf life.

PRODUCT INFO

ltem No.	Description	Content	Packaging Format
105614	Sterile DE 70% in WFI Aerosol	400ml	16 x 400ml individually double bagged
9709	Sterile IPA 70% in WFI Aerosol	400ml	16 x 400ml individually double bagged
402020	Grip Trigger	Each	1 x Grip Trigger (Non-sterile)





For validations see page 8



AGMA Healthcare Guide | 9