



## ALCOHOL 70% WFI T/SPRAYS

105665 – STERILE DE 70% WFI – 6 X 900ML T/SPRAY

0028 – NON-STERILE DE 70% WFI – 14 X 1 L T/SPRAY

1139 – NON-STERILE DE 70% WFI – 20 X 500ML T/SPRAY

DE 70% WFI IS A PROVEN MEANS FOR CONTINUAL DISINFECTION WITHIN PHARMACEUTICAL & HOSPITAL CLEANROOMS.

THE STERILE BAG-IN-BOTTLE TRIGGER SPRAYS ARE COMMONLY USED IN CLEANROOMS BOTH IN HOSPITALS AND PHARMACEUTICAL INDUSTRY THAT ARE A HIGHLY COST-EFFECTIVE MEANS OF ROUTINE DISINFECTION WHEN USED WITH STERILE DRY WIPES.

AS WELL AS A PROVEN DISINFECTANT DE 70% EVAPORATES FULLY SO THERE IS NO RESIDUE.

- EFFICACY INDEPENDENTLY VALIDATED.
- WFI QUALITY WATER – HIGH PURITY WATER.
- ERGONOMIC TRIGGER WITH EFFECTIVE RETURN MECHANISM.
- RELIABLE BOTTLE WITH HEADROOM FOR EXPANSION.
- BAG-IN-BOTTLE TO MAINTAIN STERILITY IN USE (STERILE VERSION).
- 13 WEEK IN-USE VALIDATION (STERILE VERSION).
- STERILE PRODUCT MANUFACTURED IN CLASS 100 CLEANROOMS.
- NON-STERILE PRODUCT MANUFACTURED IN CLASS 10,000 CLEANROOMS.
- FULL TRACEABILITY WITH SUPPORTING DOCUMENTATION.
- GAMMA IRRADIATED FOR STERILITY.
- DOUBLE BAGGED (STERILE VERSION).
- CONSISTENT SPRAY PATTERN



Use Biocides Safely. Always read the label and product information.

THE INFORMATION CONTAINED IN THIS DATA SHEET IS ACCURATE TO THE BEST OF OUR KNOWLEDGE AND BELIEF. NO WARRANTY IS, HOWEVER, GIVEN OR IMPLIED. CONSULT AGMA LTD CONCERNING ALL NON-STANDARD APPLICATIONS.

ISSUE NO. 2-10/22