



# STERILE ALKALINE DETERGENT

## ITEM NOS

0092 STERILE ALKALINE DETERGENT (2 X 5L)

0094 STERILE ALKALINE DETERGENT (6 X 900ML)

0096 STERILE ALKALINE DETERGENT (15 X 1L)

FOR EFFECTIVE DISINFECTION TO BE ACHIEVED IT IS IMPERATIVE THAT SURFACES ARE FREE FROM ANY VISIBLE SOILING OR BUILD-UP OF RESIDUE FROM PREVIOUS DISINFECTANT APPLICATIONS.

AGMA HAVE CONSIDERABLE EXPERIENCE IN SURFACTANT TECHNOLOGY AND CLEANING / DISINFECTION. WE HAVE COMBINED THIS EXPERTISE TO DEVELOP A STERILE ALKALINE DETERGENT SPECIFICALLY FOR USE IN CLEANROOMS AND PHARMACIES.

WITH HIGHLY EFFECTIVE SURFACE CLEANING PROPERTIES IT IS IDEAL WHEN A DEEP CLEAN IS REQUIRED. FORMULATED TO BE LOW FOAMING IT DOES NOT LEAVE SIGNIFICANT RESIDUES OR SMEARS.

- HIGHLY EFFECTIVE SURFACE CLEANING PROPERTIES.
- REMOVES GREASES, SPILLS AND DISINFECTANT RESIDUE.
- LOW FOAMING.
- DOES NOT LEAVE SIGNIFICANT RESIDUES.
- MADE WITH WFI QUALITY WATER.
- MANUFACTURED IN CLASS 10,000 CLEANROOMS.
- QUALITY MANAGEMENT SYSTEM CERTIFIED TO BS EN ISO 9001: 2015
- MANUFACTURED IN ACCORDANCE WITH cGMP GUIDELINES.
- FULLY TRACEABLE WITH SUPPORTING ANALYSIS.
- AVAILABLE IN 1L AND 5L CAP BOTTLES, AND 900ML TRIGGER BOTTLE



## DIRECTION FOR USE

AGMA STERILE ALKALINE DETERGENT IS MANUFACTURED FOR USE ON CLEANROOM HARD SURFACES TO REMOVE RESIDUES. AN INITIAL DEEP CLEAN SHOULD BE PERFORMED PRIOR TO THE USE OF ANY DISINFECTANT. THIS WILL HELP THE REMOVAL OF ANY SOILS OR HEAVY RESIDUE BUILD UP. STERILE ALKALINE DETERGENT IS AN IDEAL PRODUCT THAT CAN BE USED BETWEEN ROTATING OF DISINFECTANTS.

THE INFORMATION CONTAINED IN THIS DATA SHEET IS ACCURATE TO THE BEST OF OUR KNOWLEDGE AND BELIEF. NO WARRANTY IS, HOWEVER, GIVEN OR IMPLIED. CONSULT AGMA LTD CONCERNING ALL NON-STANDARD APPLICATIONS.

ISSUE NO. 1-05/24