



4YPIN STERILE NEUTRAL DETERGENT (DB) IN WFI PRE-IMPREGNATED – READY TO USE ITEM 0130 (10 X 50 WIPES)

A PROVEN MEANS FOR WIPING INTO PHARMACEUTICAL & HOSPITAL CLEANROOMS.
FOR EFFECTIVE DISINFECTION TO BE ACHIEVED IT IS IMPERATIVE THAT SURFACES ARE FREE FROM ANY VISIBLE SOILING OR BUILD-UP OF RESIDUE FROM PREVIOUS DISINFECTANT APPLICATIONS.

AGMA HAVE CONSIDERABLE EXPERIENCE IN SURFACTANT TECHNOLOGY AND STERILE CLEANING / DISINFECTION. WE HAVE COMBINED THIS EXPERTISE TO DEVELOP A STERILE NEUTRAL DETERGENT SPECIFICALLY FOR USE IN CLEANROOMS AND PHARMACIES.

- HIGHLY EFFECTIVE SURFACE CLEANING PROPERTIES.
- REMOVES GREASES, SPILLS AND DISINFECTANT RESIDUE.
- LOW FOAMING.
- DOES NOT LEAVE SIGNIFICANT RESIDUES.
- FORMULATED TO BE PH NEUTRAL BEFORE IRRADIATION.
- MANUFACTURED IN GRADE C/D CLEANROOMS.
- MADE WITH WFI QUALITY WATER.
- POLYPROPYLENE Z FOLDED, 30 GSM - 25CM X 25CM WIPES
- 50 WIPES PER PACK, 10 PACKS PER CASE
- GAMMA IRRADIATION FOR STERILITY.
- 2 YEAR UNOPENED SHELF LIFE.
- DOUBLE BAGGED.
- FULLY TRACEABLE WITH SUPPORTING ANALYSIS.



DIRECTION FOR USE

AGMA NEUTRAL DETERGENT WIPES ARE MANUFACTURED TO HELP CLEANROOM OPERATORS MEET THE STANDARDS OF ANNEX 1. IDEAL FOR USE ON CLEANROOM HARD SURFACES TO REMOVE RESIDUES AND FOR USE BETWEEN ROTATING OF DISINFECTANTS. AN INITIAL DEEP CLEAN SHOULD BE PERFORMED PRIOR TO THE USE OF ANY DISINFECTANT TO HELP THE REMOVAL OF ANY SOILS OR HEAVY RESIDUE BUILD UP.

THE INFORMATION CONTAINED IN THIS DATA SHEET IS ACCURATE TO THE BEST OF OUR KNOWLEDGE AND BELIEF. NO WARRANTY IS, HOWEVER, GIVEN OR IMPLIED. CONSULT AGMA LTD CONCERNING ALL NON-STANDARD APPLICATIONS.

ISSUE NO. 1-01/26