



# STERILE NEUTRAL DETERGENT MOP WIPES

## IN WFI PRE-IMPREGNATED – READY TO USE

### ITEM NO. 0085

FOR EFFECTIVE DISINFECTION TO BE ACHIEVED IT IS IMPERATIVE THAT SURFACES ARE FREE FROM ANY VISIBLE SOILING OR BUILD-UP OF RESIDUE FROM PREVIOUS DISINFECTANTS.

READY TO USE STERILE NEUTRAL DETERGENT MOP WIPES ARE HIGHLY EFFECTIVE AT REMOVING GREASES, SPILLS AND DISINFECTANT RESIDUE FROM LARGER SURFACE AREAS.

OUR PRESATURATED SINGLE-USE MOP WIPES ARE READY TO USE AND PERFECT FOR A BUCKETLESS SYSTEM, SAVING OPERATOR TIME. THE MOP WIPES ARE EASY TO INSTALL AND DISPOSE OF AND ARE COMPATIBLE WITH WIDELY AVAILABLE 'PUSH CLIP' MOP HEADS.

- REMOVES GREASES, SPILLS AND DISINFECTANT RESIDUE.
- LOW FOAMING.
- NO SIGNIFICANT RESIDUES.
- FORMULATED TO BE PH NEUTRAL BEFORE IRRADIATION.
- MADE WITH WFI QUALITY WATER.
- POLYCELLULOSE 60G HIGH QUALITY, LOW LINTING MOP WIPE
- 25 X 50CM, ALLOWING FOR LARGE AREA COVERAGE.
- FULLY TRACEABLE WITH SUPPORTING ANALYSIS.
- GAMMA IRRADIATION FOR STERILITY.
- 20 MOP WIPES PER PACK, 5 PACKS PER CASE.
- TWO YEAR UNOPENED SHELF LIFE.



Use Biocides Safely. Always read the label and product information.

THE INFORMATION CONTAINED IN THIS DATA SHEET IS ACCURATE TO THE BEST OF OUR KNOWLEDGE AND BELIEF. NO WARRANTY IS, HOWEVER, GIVEN OR IMPLIED. CONSULT AGMA LTD CONCERNING ALL NON-STANDARD APPLICATIONS.

ISSUE NO. 1-11/24